





The LexBal Intra gastric Balloon is used to reduce overweight in persons till 30 kg .Is installed by endoscopy technique, the balloon it will be filled with saline solution with methylene blue and give to the user the feel of full nest with a little amount of food the installation procedure time takes 20 minutes and the balloon it is design to be used for 6 months in the stomach ones the time is gone the extraction procedure it takes another 20 minutes in both cases the patient go to their home after one hour and continuo the live stile after 48hs. The balloon is built in implantable silicone.

The Lexel Gastric Band it is a surgical procedure made by laparoscopic technique for persons who their BMI (Body mass indicator) is over 35. It is a very simple procedure well know in the morbid obesity market. One of the big advantages over other procedures are that the banding is reversible procedure that mean that when the patient have their new alimentary habits and have the right weight the device could be take it out. It is Built in Implantable silicone with and implanted port to where is manage the flow by the stomach.

The TITANIUM-PORT implanted port it is used for long term chemotherapy. It could be implanted till 5 years and resist around 1500 punctures with a 22 G Huber needle. The port it is built in Titanium and silicone those materials allow the use of MRI and CTA. The set comes with a complete line of accessories to make easy their implantation by the different technique of implantation. Titanium-Port have several models Adult, Pediatric and Neonatal and also different catheter 5,5 FR, 7 FR and 9 FR.

The POLY-PORT Huber needles are design to access an implantable Port to chemotherapy. The set comes with a clamp and side port to infusion. The Poly-Port has an anatomical design to make easier the fixation of the needles to the patient in long infusions. The Poly port have several models to all the port access requirement.

Model	Length		
PP 19G	17mm	25mm	
PP20G	17mm	25mm	32mm
PP21G	17mm	25mm	32mm
PP22G	17mm	25mm	32mm

Lexel Bone Marrow Biopsy and Aspiration needle have a hardness and Sharpe who allow an easy insertion in the Bone Marrow to take samples and make a diagnostics. We also have it in different sizes.

Model	Length	
Bone Marrow Biopsy 11G	10cm	15cm
Bone Marrow Biopsy 8G	10cm	15cm
Bone Marrow Aspiration 15G	4cm	
Bone Marrow Aspiration 16G	4cm	